

BAREFOOT FREEDOM®

2016



NEW FOR FALL 2016

HOPE

SEE PAGE 9



ASHTON

SEE PAGE 7

KOOL

SEE PAGE 9



GENOA

SEE PAGE 6

BERLIN

SEE PAGE 6





SAVANNAH SEE PAGE 9



SEE PAGE 5
EXCEL



TUSCANY SEE PAGE 4

SEE PAGE 5
SWIFT



MORGAN SEE PAGE 4



SEE PAGE 19
MEN'S DELAWARE

Athleisure



NEW! Morgan

14281-19 Black Mesh/ Stretch

HCPCS Code A5500 Morgan & Tucson pending

PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Flared, Lightweight Rubber Outsole with EVA Wedge Midsole
- Removable, Polyurethane Footbed with Heel Cushion and Maximum Built-up Arch Support
- Full Contact Rocker Bottom
- Built-In Wide, Synthetic Fiber Shank
- Firm Heel Counter
- Accommodates Prescribed Orthotics



NEW! Tuscany

10555-19 Black Suede/ Black Mesh

10555-8G Pewter Calf/ Pewter Mesh

N (AA) 7-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13
XW (4E) 5-12, 13

239



Delite

14373-19 Black Calf/Black Mesh

14373-21 White Calf/White Mesh

14373-23 Sky Blue Nubuck/White Mesh

14373-31 Bone Nubuck/Bone Mesh

14373-58 Grey Nubuck/Grey Stretch (color shown)

14373-99 Black Calf/Black Stretch

HCPCS Code A5500 Joy, Bliss, Delite

AD RI SS

Barefoot Freedom®

by DREW

- Added Depth with Removable Footbed
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control (Delite)
- Soft Fabric Lining (Joy & Bliss)
- Flexible, Lightweight, Polyurethane Outsole
- Firm Heel Counter
- Tempered Steel Shank
- Accommodates Prescribed Orthotics

N (AA) 6-12
M (B) 5-12
W (D) 5-12
WW (EE) 5-12

* Delite Not Available in 11.5

200



Bliss

10389-53 Black Suede

10389-4V Navy Suede (color shown)

Adjustable Backstrap!



Joy

14368-19 Black Smooth Leather/Black Mesh with Floral Trim

14368-21 White Smooth Leather/White Mesh with Floral Trim

14368-31 Bone Smooth Suede/ Bone Mesh with Floral Trim (color shown)

14368-41 Navy Suede/ Navy Stretch

14368-4C Sky Blue Suede/ White Mesh with Floral Trim

14368-69 Brown Suede/ Brown Stretch

14368-7C Olive Suede/ Olive Stretch

14368-99 Black Smooth Leather/ Black Stretch



NEW! Excel

10164-19 Black Sport Mesh/ Stretch
10164-41 Navy Blue Sport Mesh/ Stretch



NEW! Swift

10077-10 Black Sport Mesh
10077-90 Royal Blue Sport Mesh

HCPCS Code A5500 Solo (Excel & Swift pending);
excludes Blast



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable, Polyurethane, Drilex® Lined Insole
- Breathable Mesh and Soft Fabric Lining
- Durable Rubber Outsole and EVA Midsole
- Dual Shank System
- Wide, Synthetic Fiber Shank
- External TPU Shank
- Firm Heel Counter
- Rocker Bottom
- Accommodates Prescribed Orthotics

N (AA) 7-12*
M (B) 5-13
W (D) 5-13
WW (EE) 5-13
XW (4E) 5-13

Narrow width not Available in Blast



Solo

14214-10 Black Sport Mesh
14214-20 White Sport Mesh
14214-4F Royal Blue Sport Mesh
14214-57 Red Sport Mesh



Blast

10065-10 Black Sport Mesh

Women's

Casuals



NEW! Genoa

14316-19 Black Leather Combo
14316-41 Navy Leather Combo



NEW! Berlin

13243-12 Black Leather
13243-42 Navy Leather



NEW! London II

13252-41 Navy Leather Combo (color shown)
13252-99 Black Leather Combo



Prague

10234-12 Black Leather
10234-86 Pewter Leather (color shown)



Cairo

14491-12 Black Leather
14491-64 Brown Leather (color shown)

HCPCS Code A5500 Cairo, Prague (Genoa pending);
excludes London II and Berlin



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable Polyurethane Footbed
- Flexible, Lightweight, Polyurethane Outsole
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control (Cairo & Prague)
- Leather Linings (Genoa & Berlin)
- Soft, Fabric Lining (London)
- Drilex® lined footbed with AEGIS Microbe Shield™ and Metatarsal Support
- Tempered Steel Shank
- Firm Heel Counter
- Accommodates Prescribed Orthotics

N (AA) 6-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13

235

NEW!

Medial Zipper Closure



NEW! Ashton

13340-12 Black Leather

13340-68 Dark Brown Leather

Velcro® Closure



Summer

14420-12 Black Leather

14420-42 Navy Leather (color shown)

14420-86 Pewter Leather (color shown)



Velcro® Closure

Monaco

14435-12 Dusty Black Leather

HCPCS Code A5500 Summer & Monaco (Ashton pending)

AD RI SS

- Added Depth with Removable Footbed
- Leather Lining (Summer)
- Lightweight, Slip-resistant, Rubber Outsole with EVA Midsole
- Firm Heel Counter
- Tempered Steel Shank
- 14/8 Heel Height
- Accommodates Prescribed Orthotics

N (AA) *6-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13

Sizes 6 & 6.5 Narrow not Available in Ashton



236

HCPCS Code A5500 Jada

AD RI SS

- Added Depth
- Removable, Moldable, Dual-density Insole with Shock Absorption
- Steel Shank for Added Support
- Polyurethane Outsole
- Firm Heel Counter
- 6/8 Heel Height
- Accommodates Prescribed Orthotics

26

M (B) 6-11, 12
W (D) 6-11, 12
WW (EE) 6-11, 12




Jada

14320-3Z Ivory Stretch

14320-99 Black Stretch

Women's

Booties & Shooties & Hikers, Oh My!



Waterproof!

Glacier

10188-09 Black/Grey Nubuck with Yellow

10188-6K Brown Nubuck with Orange

HCPCS Code A5500 Glacier

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Polyurethane Footbed with Added Heel Cushion
- Waterproof Upper
- Waterproof Lining
- Slip Resistant, Oil Resistant and Non-Marking Outsole
- Durable, Lightweight Rubber Outsole
- Firm Heel Counter
- Tempered Steel Shank
- Accommodates Prescribed Orthotics

N (AA) 6-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13

495



Water Resistant!

Now Available in Smooth Leather!



Krista

10141-12 Black Leather New Color!

10141-53 Black Suede

10141-6R Brown Suede



Water Resistant!

Haley

13344-53 Black Suede

13344-6R Brown Suede

HCPCS Code A5500 Krista, Haley

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Polyurethane Footbed with Drilex® Topcover
- Flexible, Lightweight, Polyurethane Outsole
- Drilex® Lining with AEGIS Microbe Shield™
- Water Resistant
- Firm Heel Counter
- Tempered Steel Shank
- Accommodates Prescribed Orthotics

N (AA) 6-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13

234

Insoles



**Drew® ULTRON
Footbed**

Barefoot Freedom
By Drew

Comfort Insoles

- Unique micro-cellular polyurethane formulation delivers unsurpassed comfort and performance
- Contoured arch for extra support
- Stabilizing heel cradle for shock absorption
- Moisture wicking Drilex® topcover keeps your feet cool and dry
- Energy return system for maximum performance

80100 Barefoot Freedom by Drew Footbed

Fits Women's Shoe Size N M W Widths	Insole Size
5 to 6½.....	1
7 to 8½.....	2
9 to 10½.....	3
11 to 12.....	4
Fits Women's Shoe Size WW XW Widths	Insole Size
5 to 6½.....	5
7 to 8½.....	6
9 to 10½.....	7
11 to 12.....	8

ALL NEW!

Adjustable Swivel Strap

NEW! Savannah

Stretch Vamp



17100-9T Grey Flannel
17100-99 Black Stretch

NEW! Hope



10573-53 Black Suede
10573-6R Brown Suede

NEW! Kool



19178-19 Black Calf/ Black Suede

HCPCS Code A5500 Savannah, Hope & Kool Pending

(DD) (PF) (RI) (SS)

- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Polyurethane Footbed with Drilex® Topcover
- Flexible, Lightweight, Polyurethane Outsole
- Drilex® Lining with AEGIS Microbe Shield™
- Firm Heel Counter
- Tempered Steel Shank
- Accommodates Prescribed Orthotics

N (AA) 7-12
M (B) 5-12, 13
W (D) 5-12, 13
WW (EE) 5-12, 13




234

Women's

Footwear and Last Recommendations for Common Foot Problems


Common Foot Problems



Bunions: Bunions require footwear with a wide and deep toe box to accommodate the affected joint. Drew brand footwear has added depth or double depth, large forefoot capacity and soft yielding uppers. Lasts that work well for bunions include women's lasts 88, 46, 221, 237 and men's lasts 66 and 420.



Heel Pain: Heel pain – especially when standing – is helped by customized cushion inserts and added depth footwear. Heel pain is usually caused as a result of Plantar Fasciitis. Drew footwear with prescribed orthotics works well for this condition. Lasts that work well for heel pain include women's lasts 239, 801 and men's lasts 428 and 741.



Pes Planus: Correctly designed shoes will make the low-arched foot more comfortable. Depending on mobility, women's lasts 46, 88 and men's last 66, 420 and 401 work well for the rigid foot, while last 101, 401 and 420 works well for the flexible flat foot.



Plantar Fasciitis: Pain in the arch area is a symptom of Plantar Fasciitis, which can be managed by customized orthotics and proper footwear. All Drew brand footwear can accommodate prescribed orthotics, which help relieve the pain of Plantar Fasciitis. Lasts that work well for Plantar Fasciitis include women's lasts 235, 234 and men's lasts 420.



Hammertoes: Shoes with wide toe boxes and metatarsal lifts help the pain of hammertoes. Drew brand footwear with high toe boxes, removable cushion inserts, and soft leather uppers benefit this condition. Look for women's lasts 88 and 46 and men's last 66.



Ankle Valgus: Ankle roll in requires strong shoes with wedges or flares that reduce pain in the lower legs and ankles. Drew brand footwear, with its supportive heel counter, wide shank area, rigid heel tuck, and flared outsoles are ideal for Ankle Valgus. This footwear offers complete rearfoot control, thus helping balance and align the foot. Lasts that work well for Ankle Valgus include women's lasts 239, 101 and men's lasts 66 and 401.



Calluses: Accommodative inserts to reduce pressure will allow most calluses to lessen and sometimes disappear. Drew brand footwear, constructed with moldable, cushion inserts are accommodative to this condition.



Pes Cavus: A high-arched foot requires therapeutic footwear. A long medial counter provides support and stability. Drew lasts with modified inserts or orthotics are also beneficial for this problem. Lasts that work well for Pes Cavus include women's last 235 and men's last 66.



Ankle Varus: Ankle roll out requires strong footwear with wedges or flares that reduce pain in the lower legs and ankles. Drew footwear with its supportive counter, wide shank areas, rigid heel tuck, and flared outsoles are ideal for Ankle Varus. This footwear offers complete rearfoot control, thus helping balance and align the foot. Lasts that work well for Ankle Varus include women's lasts 239, 101 and men's lasts 66 and 420.

Last Recommendations for Common Foot Problems

Women's

Pedic 221 – See Pages 12 and 13

- Added depth combined with broad forepart
- Wide shank area for superior lateral stability
- Excellent for application to treat minor foot conditions
- 6/8 heel height

231 & 237

Pedic 200 – See Page 4

- Snug fitting heel seat
- Short back last
- Excellent for applications to treat minor foot conditions
- 6/8 heel height

200

Pedic 234 – See Pages 4, 8 and 9

- Added depth combined with broad forepart
- Wide base with mid-foot stabilizing shank provide arch support and stability
- Excellent last to accommodate prescribed orthotics
- 8/8 heel height

234

Pedic 80 & 801 – See Pages 5 and 14

- Added depth combined with broad forepart
- Wide shank provides superior lateral stability
- Excellent last to accommodate prescribed orthotics

80 & 801

Pedic 231 – See Page 11

- Added depth combined with broad forepart
- Wide base with mid-foot stabilizing shank provide arch support and stability
- Excellent last to accommodate prescribed orthotics
- 8/8 Heel Height

231

Pedic 101 – See Page 15

- Wide shank area provides superior lateral stability and also makes this last extremely versatile for orthotic applications
- Deep, broad forepart offers maximum capacity

101

Pedic 16 – See Page 17

- Added depth combined with broad forepart
- Wide shank area for superior lateral stability
- Excellent for application to treat minor foot conditions
- 6/8 heel height

16

Pedic 88 – See Page 16

- A universal last used in many applications for extreme foot conditions
- Deep, broad, oblique forepart offers more capacity than any of our other women's lasts
- Wide shank area provides superior lateral stability – extremely versatile for orthotic applications
- 8/8 heel height

88

Pedic 46 – See Page 16

- Wide shank area provides superior lateral stability and also makes this last extremely versatile for orthotic applications
- Deep, broad, oblique forepart offers maximum capacity
- 6/8 heel height offers excellent balance and weight distribution

46

Pedic 495 – See Pages 8 & 24 (Women's & Men's)

- Added depth combined with broad forepart
- Wide shank area for superior lateral stability
- Excellent for application to treat minor foot conditions
- 6/8 heel height

495

Pedic 235 – See Page 6

- Added Depth
- Universal Last
- Excellent for applications to treat minor foot conditions
- 8/8 heel height

235

Men's

Pedic 420 – See Page 24

- Deep, broad, oblique toe offers maximum capacity
- Wide shank area provides superior lateral stability
- Excellent last to accommodate prescribed orthotics

420

Pedic 66 – See Pages 20 & 21

- Snug fitting heel with extra displacement in cuboid
- Wide shank area provides excellent lateral stability
- Excellent last to accommodate prescribed orthotics
- Additional overall capacity

66

Pedic 401 – See Page 22

- Wide shank area provides superior lateral stability and also makes this last extremely versatile for orthotic applications
- Deep, broad forepart offers maximum capacity
- 8/8 heel height

401

Pedic 74 & 741 – See Page 23

- Added depth combined with broad forepart
- Wide shank provides superior lateral stability
- Excellent last for various orthotic applications
- 8/8 heel height

741

Pedic 428 – See Page 19

- Added depth combined with broad forepart
- Wide shank provides superior lateral stability
- Excellent last to accommodate prescribed orthotics

428

Stretch Vamp
& Arch



Stretch Vamp
& Arch*



Amora

13380-98	Brown Nubuck Combo
13380-99	Black Nubuck Combo

Antwerp

14331-15	Black Nubuck
14331-3H	Stone Nubuck (color shown)
14331-65	Brown Nubuck
14331-77*	Taupe Leather/ Taupe Stretch
14331-99*	Black Leather/ Black Stretch (color shown)



Geneva

10589-1Q	Black Leather/Nubuck
10589-43	Grey Leather/ Nubuck (color shown)
10589-4Q	Blue Denim Leather/Nubuck
10589-6Z	Espresso Leather/Nubuck



Lisbon

10679-19	Black Nubuck Combo
10679-69	Brown Nubuck Combo

HCPCS Code A5500 Geneva

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable ULTRON® Footbed for EVERLASTING Comfort
- Moldable Rubber Outsole with Cork
- Drilex® Upper and Footbed Lining with AEGIS Microbe Shield
- Soft, Full Grain Leather Uppers
- Tempered Steel Shank
- Extended Medial Heel Stabilizer
- Accommodates Prescribed Orthotics

223

N (AA)	6-12
M (B)	5-12
W (D)	5-12
WW (EE)	5-12
XW (4E)	5-11

HCPCS Code A5500 Antwerp, Lisbon (excludes Amora)

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable Polyurethane Footbed with Drilex® Topcover and Added Heel Cushion and Metatarsal
- Ultra Lightweight Polyurethane Outsole
- Extended Medial Heel Stabilizer
- Tempered Steel Shank
- Soft, Fabric Lining
- Accommodates Prescribed Orthotics

231

N (AA)	6-12
M (B)	5-12, 13
W (D)	5-12, 13
WW (EE)	5-12, 13
XW (4E)	5-12, 13

Drew Shoe uses genuine Velcro®. Velcro® is a registered trademark of Velcro Industries B.V.

Women's

Active



Tulip

10202-12	Black Full Grain Leather^*
10202-1P	Black Croc Patent Leather^*
10202-22	White Full Grain Leather^
10202-32	Bone Pebble Leather
10202-64	Brown Full Grain Leather
10202-6P	Brown Croc Patent Leather
10202-9Q	Copper Metallic Leather^ (color shown)
10202-99	Black Calf/ Black Stretch*#



Magnolia

14326-41	Navy Nappa/ Stretch*
14326-69	Brown Nappa/ Stretch* (color shown)
14326-99	Black Nappa/ Stretch*

HCPCS Code A5500 Tulip, Rose, Magnolia, Lotus, Orchid, Ginger, Heather



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable ULTRON® footbed with Drilex® Topcover
- Ultra Lightweight Polyurethane Slip-resistant Outsole
- Extended Medial Heel Stabilizer
- Tempered Steel Shank
- Drilex® Lining with AEGIS Microbe Shield
- Accommodates Prescribed Orthotics

^S (4A) 6-12
N (AA) 6-13#
M (B) 5-13*
W (D) 5-13*
WW (EE) 5-13*

(^Indicates Available Colors in Slim Widths^)

(*Indicates Available in sizes 12.5 & 13 Medium, Wide and Wide Wide*)

(#Indicates Available in sizes 12.5 & 13 Narrow#)

(*Indicates Available in sizes 13 Medium, Wide and Wide Wide*)



221



Rose

14375-1P	Black Croc Patent Leather^*
14375-42	Navy Leather (color shown)
14375-48	White Floral Snake Print Leather
14375-57	Red Full Grain Leather
14375-69	Brown Full Grain Leather
14375-6P	Brown Croc Patent Leather*
14375-86	Dusty Pewter Leather (color shown)
14375-8F	Brown Print Patent Leather
14375-8Y	Pewter Metallic Leather
14375-99	Black Full Grain Leather^*
14375-9L	Bone Pebble Leather^*
14375-9Q	Copper Metallic Leather^



Lotus

14631-12	Black Full Grain Leather*
14631-22	White Full Grain Leather
14631-32	Bone Full Grain Leather (color shown)



Orchid

14370-12	Black Full Grain Leather*
----------	---------------------------

Active



237

Ginger

- 14321-12 Dusty Black Leather`
- 14321-86 Dusty Pewter Leather`
- 14321-9G Dusty Gold Leather`



237

Heather

- 14400-12 Black Leather`
- 14400-83 Brandy Leather`

Women's

Athletic



Athena

10268-02	White Leather/Mesh with Grey
10268-04	White Leather/Mesh with Blue
10268-09	Black Leather/Mesh with Silver
10268-3V	Cream Suede/Mesh (color shown)

HCPCS Code A5500 Athena

DD PF RI



80

- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Flared, Lightweight Rubber Outsole with EVA Wedge Midsole
- Removable, Polyurethane Footbed with Heel Cushion and Arch Support
- Firm Heel Counter
- Wide, Synthetic Fiber Shank
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

N (AA)	7-12
M (B)	5-12
W (D)	5-12
WW (EE)	5-12
XW (4E)	5-12



Fusion

10545-12	Black Leather
10545-22	White Leather

HCPCS Code A5500 Fusion, Flash II, Flash II V, Paige

DD PF RI



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable Polyurethane, Drilex® Lined Footbed
- Breathable Mesh and Soft Fabric Lining
- Durable Rubber Outsole and EVA Midsole
- Dual Shank System
- Wide, Synthetic Fiber Shank
- External TPU Shank
- Firm Heel Counter and Rocker Bottom
- Accommodates Prescribed Orthotics

801



Flash II

10560-07	White Leather / Grey Mesh (color shown)
10560-19	Black Leather/ Black Mesh
10560-21	White Leather/ White Mesh



Flash II V

14565-07	White Leather / Grey Mesh
14565-19	Black Leather/ Black Mesh
14565-21	White Leather/ White Mesh



Paige

14695-12	Black Leather^
14695-22	White Leather^
14695-32	Bone Leather^

N (AA)	7-13*
M (B)	5-13
W (D)	5-13
WW (EE)	5-13
XW (4E)	5-13

*Narrow Size 12.5 & 13 in Paige Only

Drew Shoe uses genuine Velcro®.
Velcro® is a registered trademark of Velcro Industries B.V.

Women's

Athletic & Sandals



Motion

10186-12	Black Leather
10186-22	White Leather



Flare

10285-04	White/Blue Leather/White Mesh (color shown)
10285-21	White/White Perf Leather/White Mesh
10285-31	Bone Leather/Bone Mesh



Motion V

14406-12	Black Leather
10186-22	White Leather



Toggle Closure with Stretch Laces

Element

DD RI

17031-15	Black Nubuck
17031-4R	Denim Nubuck (color shown)
17031-6S	Brown Nubuck

233

Cascade

AD RI

17051-15	Black Nubuck
17051-2F	Sport White Leather
17051-2R	Sand Nubuck (color shown)
17051-4R	Denim Nubuck



HCPCS Code A5500 Element (excludes Cascade)

- Double Depth (Element)
- Added Depth (Cascade)
- EVA & Memory Foam Footbed with Sueded Topcover (Element)
- Polyurethane Footbed with Sueded Topcover and Metatarsal Support (Cascade)
- Soft Fabric Lining
- Lightweight, Rubber Outsole with EVA Wedge Midsole
- Tempered Steel Shank
- Velcro® Brand Fastened Footbed
- Accommodates Prescribed Orthotics

N (AA)	6-12, 13
M (B)	5-13
W (D)	5-13
WW (EE)	5-13

Removable
INSOLE

230

HCPCS Code A5500 Motion, Motion V, Flare

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Tri-support Insole with Built-In Heel Cushion Pad
- Rocker Bottom and Wide Shank
- Durable Rubber Outsole with EVA Midsole
- Tempered Steel Shank
- Extended Medial Heel Stabilizer
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

N (AA)	7-12
M (B)	5-13
W (D)	5-13
WW (EE)	5-13
XW (4E)	5-13

101

Ariana

17488-2B	White/Pink Floral Leather
17488-86	Pewter Leather



Alana

17060-1R	Black & White Leather
17060-29	Orange & Yellow Leather*
17060-8Z	Pewter & Bronze Leather
17060-9K	Green & Fuchsia Leather (color shown)



DD RI SS

- Double Depth
- Thick, Removable ULTRON Footbed with Memory Foam and Sueded Topcover
- Full Leather Linings
- Lightweight, High-abrasion, Dual-density, Durable EVA Outsole
- Tempered Steel Shank
- Velcro® Brand Fastened Footbed
- Accommodates Prescribed Orthotics

Removable
INSOLE

120

*N (AA)	6-13
M (B)	5-13^
W (D)	5-13^
WW (EE)	5-13^

*No Narrow Width Available
^Size 12.5 Available in Alana Only

Women's

Therapeutic

Blazer

10418-12 Black Calf

10418-75 Taupe Nubuck

N (AA) 7-11, 12
M (B) 6-11, 12
W (D) 6-11, 12
WW (EE) 6-11, 12
XW (4E) 6-11



New Villager

10676-11 Black Soft Pebble Leather

N (AA) 7-11, 12
M (B) 5-11, 12
W (D) 6-11, 12
WW (EE) 6-11, 12
XWW (4E) 6-10



Extra Wide
Opening for Easy
On- Easy Off

Seamless, Moldable
Deerskin Linings

Phoenix Plus

10217-12 Black Calf

S (AAA) 7-11, 12
N (AA) 7-11, 12
M (B) 5-11, 12
W (D) 6-11, 12
WW (EE) 6-11, 12

HCPSC Code A5500 Blazer, New Villager, Victoria, Phoenix Plus

DD PF RI SS PLUS

- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable, Moldable, Dual-Density Insert
- Drilex® Lining with AEGIS Microbe™ Shield for Odor and Stain Control
- Flared Outsole with Polyurethane Rocker Bottom
- Extended Medial Heel Stabilizer
- Tempered Steel Shank
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

46

Victoria

14463-12 Black Full Grain Leather

M (B) 6-11, 12
W (D) 6-11, 12
WW (EE) 6-11, 12

HCPSC Code A5500 Fiesta, Sedona

AD RI SS

- Added Depth
- Removable, Moldable, Dual-density Insert
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Lightweight Polyurethane Outsole
- Rocker Bottom
- Tempered Steel Shank
- Extended Medial Heel Stabilizer
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

N (A) 6-12
M (C) 5-12
W (E) 5-12
WW (EEE) 5-12
XWW (5E) 5-11

88

Fiesta

10175-12 Black Calf



Sedona

10125-12 Black Calf



Women's

Stretch Vamp



Quest

14251-99 Black Leather/ Black Stretch*



Bloom II

14353-12	Black Calf
14353-32	Bone Calf
14353-3L	Bone Calf/Bone Stretch
14353-42	Navy Calf
14353-75	Taupe Nubuck (color shown)
14353-99	Black Calf/Black Stretch



Parade II

10295-12	Black Calf
10295-22	White Calf
10295-32	Bone Calf
10295-72	Taupe Calf (color shown)

Paradise II

14521-12	Black Calf
14521-22	White Calf
14521-32	Bone Calf (color shown)
14521-72	Taupe Calf



HCPCS Code A5500 Quest, Parade II, Zip II, Bloom II, Paradise II

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable, Moldable, Dual-density Insert
- Flared Outsole with Wide Shank
- Polyurethane Rocker Bottom
- Firm Heel Counter
- Tempered Steel Shank
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Soft, Foam-padded Collar
- Accommodates Prescribed Orthotics

*N (AA)	6-12
M (B)	5-12
W (D)	5-12
WW (EEE)	5-12

*Quest not available in Narrow Width

16



Zip II

10821-18 Black Tumbled Nappa



Footsaver® Medical Grade Inserts

80600 Orthotic

HCPCS Code A5512

Drew's most advanced, comfort-based orthotic. Ready for quick and easy modification.

- Drew-Lite® Cushioning Topcover
- EVA Foot Flex Base 3/16" Shore A 40 Durometer

Insoles

Fits Women's Shoe Size	Fits Men's Shoe Size	Insole Size
5 to 6½	N/A	1
7 to 8½	6 to 7½	2
9 to 10½	8 to 9½	3
11 and up	10 to 11½	4
N/A	12 and up	5

Women's

Drew Shoe Solutions for Common Foot Problems

Exceptional footwear that helps alleviate common foot problems.

This quick reference guide offers helpful information for those fitting and selling Drew Shoes.

Key to Shoe Ratings for Specific Foot Problems

1

The **BEST** shoe for this particular problem

2

A **GOOD** shoe for this particular problem

3

An **ACCEPTABLE** shoe for this particular problem with the necessary orthotic or pedorthic modification

WOMEN'S

	Page Number	Last Number	Pes Cavus (High Arch)	Pes Planus (Flat Foot)	Ankle Varus (Roll Out)	Ankle Valgus (Roll In)	Metatarsalgia (Ball of Foot Pain)	Plantar Fascitis (Heel Pain)	Bunions	Hammertoes	Arthritis	Diabetes
Morgan	4	239	2	1	1	1	2	2	2	2	2	2
Tuscany	4	239	2	1	1	1	2	2	2	2	2	2
Delite	4	200	3	2	3	2	3	3	2	3	3	3
Bliss	4	200	2	2	2	2	2	2	2	2	2	2
Joy	4	200	3	2	3	2	3	2	2	2	2	2
Swift	5	801	2	2	2	2	1	2	1	1	1	2
Excel	5	801	2	2	2	2	1	2	1	1	2	2
Solo/ Blast	5	801	2	2	2	2	1	2	1	2	1	3
Genoa	6	235	2	2	2	2	2	2	3	3	2	3
Berlin	6	235	2	2	2	2	2	2	3	3	2	3
London II	6	235	2	2	2	2	2	2	2	2	2	3
Cairo	6	235	1	1	2	2	2	2	1	2	1	2
Prague	6	235	1	1	2	2	2	2	2	2	2	2
Ashton	7	236	2	2	2	2	3	2	3	2	2	2
Summer/ Monaco	7	236	2	2	2	2	3	2	2	3	2	3
Krista	8	234	2	2	2	2	2	2	1	2	2	2
Haley	8	234	2	2	2	2	2	2	1	2	1	2
Savannah	9	234	2	2	3	3	2	2	1	2	2	2
Hope	9	234	2	2	2	2	2	2	1	2	2	2
Kool	9	234	2	2	2	2	2	2	1	2	2	2
Amora/ Antwerp	11	231	2	2	2	2	2	2	1	1	1	2
Lisbon	11	231	2	2	2	2	2	2	1	1	2	2
Geneva	11	223	2	2	2	1	1	2	2	2	2	2
Tulip	12	221	2	2	2	2	2	2	1	2	2	2
Magnolia	12	221	2	2	2	2	2	2	1	1	1	1
Rose	12	221	2	2	2	2	2	2	1	2	1	2
Lotus/ Orchid	12	221	2	2	2	2	2	2	1	2	1	2
Ginger	13	221	2	2	2	2	2	2	2	2	2	2
Heather	13	221	2	2	2	2	2	2	2	2	1	2
Athena	14	80	2	2	2	2	1	2	2	2	3	3
Paige	14	801	2	2	2	2	1	2	2	2	1	3
Fusion	14	801	2	2	2	2	1	2	2	2	2	3
Flash II/Flash II V	14	801	2	2	2	2	1	2	2	2	2	3
Motion/ Motion V	15	101	2	1	1	1	1	2	3	3	3	3
Flare	15	101	2	1	1	1	1	2	3	3	3	3
Blazer	16	46	3	1	1	1	1	2	2	2	2	2
New Villager	16	46	3	1	1	2	1	2	1	1	1	1
Victoria	16	46	3	1	1	2	1	2	1	1	1	2
Phoenix Plus	16	46	3	1	1	1	1	2	2	2	2	2
Fiesta	16	88	2	1	1	1	1	2	1	1	1	1
Sedona	16	88	2	1	1	1	1	2	1	1	1	1
Quest	17	16	2	2	2	2	2	3	1	1	1	1
Bloom II	17	16	3	2	3	3	2	3	3	3	2	3
Parade II	17	16	2	2	2	2	2	3	2	2	2	3
Paradise II	17	16	2	2	2	2	2	3	2	2	2	3
Zip II	17	16	2	2	2	2	2	3	2	2	2	2

MEN'S

Mansfield	19	428	2	2	2	2	1	2	2	2	1	2
Bexley/Toledo	19	428	2	2	2	2	1	2	2	2	2	2
Dakota	19	428	2	2	2	2	1	2	2	2	1	2
Keith	19	428	2	2	2	2	1	2	2	2	2	2
Delaware	19	428	2	2	2	2	1	2	2	2	2	2
Navigator II	20	66	2	2	2	2	1	2	2	2	2	3
Expedition II	20	66	2	2	2	2	1	2	2	2	2	3
Big Easy/ Tucson	20	66	2	2	2	2	1	2	1	1	2	2
Traveler/Traveler V	21	66	2	1	1	1	1	2	2	2	2	2
Victor	21	66	3	1	1	2	1	2	1	1	1	2
Walker II	21	66	2	2	2	2	1	2	2	2	2	3
Journey II	21	66	2	2	2	2	2	3	1	1	1	1
Doubler	21	66	2	1	1	1	1	2	1	1	1	1
Force/ Force V	22	401	2	1	1	1	1	2	3	3	3	3
Pulse/ Surge	22	401	2	1	1	1	1	2	3	3	3	3
Voyager	18	741	2	2	2	2	1	2	2	2	2	3
Lightning II/ V	18	741	2	2	2	2	1	2	2	2	2	3
Jimmy	18	741	2	2	2	2	1	2	2	2	2	3
Aaron	19	74	2	2	2	2	1	2	2	2	3	3
Trail/ Trek	23	495	2	2	2	2	2	2	2	3	3	3
Rockford	23	420	2	1	1	1	2	2	1	2	2	2
Boulder	23	420	2	1	1	1	2	2	1	2	2	2

BAREFOOT FREEDOM®

To Order:
phone 02- 94385554 fax
02- 94385571
1-657 Pacific Highway St.
Leonards NSW 2065
www.barefootfreedom.com.au

Men's Casuals

**NEW
Color!**



Toledo

40895-14	Black Leather
40895-65	Brown Nubuck (color shown) New Color!
40895-68	Brown Leather
40895-83	Brandy Leather

NEW!



NEW! Delaware

40939-6R	Dark Brown Suede
40939-9T	Dark Grey Suede



Mansfield

44915-14	Black Leather
44915-68	Brown Leather



**Twin Gore
Stretch**

Bexley

43919-14	Black Leather
43919-60	Brown Leather (color shown)
43919-65	Brown Nubuck
43919-80	Tan Leather



Keith

40672-14	Black Leather
40672-83	Brandy Leather

HPCPS Code A5500 Mansfield, Dakota, Toledo, Keith (excludes Bexley); Delaware Pending

DD PF RI SS PLUS

- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable Polyurethane Footbed
- Opanka Construction
- Drilex® and Leather Lining
- Tempered Steel Shank
- Lightweight Polyurethane Outsole with Rocker Bottom
- Firm Heel Counter
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

428



Dakota

40475-14	Black Leather
40475-83	Brandy Leather

N (B)	9-12, 13
M (D)	8-12, 13, 14, 15, 16
W (EE)	8-12, 13, 14, 15, 16
4W (4E)	7-12, 13, 14, 15, 16
6E	7-12, 13, 14, 15, 16

Drew Shoe uses genuine Velcro®.

To Order: phone 02- 94385554 Fax 02- 94385571 www.barefootfreedom.com.au

Men's

Men's Sandal & Therapeutic



Dublin

47717-2U Olive Suede

499

M (D) 8-13, 14, 15, 16
W (EE) 8-13, 14, 15, 16
4W (4E) 7-13, 14, 15, 16
6E 7-13, 14, 15, 16



Warren

47791-09 Black / Grey Nubuck

47791-84 Brown / Tan Nubuck



Mason

47726-18 Black Tumbled Leather

47726-65 Brown Nubuck

M (D) 8-13, 14, 15, 16
W (EE) 8-13, 14, 15, 16
4W (4E) 7-13, 14, 15, 16
6E 7-13, 14, 15, 16

HCPSC Code A5500 Dublin pending; excludes Warren & Mason

DD RI SS

- Double Depth
- Removable, EVA Footbed with Memory Foam (Dublin)
- Removable, Cushioned ULTRON Footbed (Warren & Mason)
- Sueded Lining (Dublin)
- Soft Stretchable Lining (Warren & Mason)
- Flexible, Shock-absorbing, Rubber Outsole
- Tempered Steel Shank
- Velcro® Brand Fastened Footbed
- Accommodates Prescribed Orthotics
- Adjustable Velcro® Brand Closures



Navigator II

44837-14 Black Pebbled Leather

44837-64 Brown Pebbled Leather

N (B) 8½-12, 13, 14
M (D) 8-12, 13, 14, 15, 16
W (EE) 7-12, 13, 14, 15, 16
4W (4E) 7-12, 13, 14, 15, 16
6E 7-12, 13, 14, 15, 16



Big Easy

44859-14 Black Pebbled Leather

44859-68 Brown Pebbled Leather



Tucson

40678-14 Black Pebbled Leather

N (B) 9-13, 14
M (D) 8-13, 14, 15, 16
W (EE) 8-13, 14, 15, 16
4W (4E) 7-13, 14, 15, 16
6E 7-12, 13, 14, 15, 16

HCPSC Code A5500 Navigator II, Big Easy, Tucson

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Moldable, Dual-density, Contoured, Polyurethane Cushion Footbed
- Extended Medial Heel Stabilizer
- Lightweight, Flared Polyurethane Outsole
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Tempered Steel Shank
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

Men's

Men's Therapeutic



Traveler

40973-14 **Black Calf**
 40973-59 **Cognac Nubuck**
 N (B) 9-12, 13
 M (D) 8-12, 13, 14, 15
 W (EE) 8-12, 13, 14, 15
 4W (4E) 7-12, 13, 14
 6E 7-12, 13



Traveler V

44908-14 **Black Calf**
 44908-59 **Cognac Nubuck**
 M (D) 8-12, 13, 14, 15
 W (EE) 8-12, 13, 14, 15
 4W (4E) 7-12, 13, 14
 6E 7-12, 13



Victor

44700-14 **Black Pebbled Leather**
 M (D) 8-12, 13, 14, 15
 W (EE) 8-12, 13, 14, 15
 4W (4E) 8-12, 13, 14



Journey II

44885-99 **Black Leather/Black Stretch**
 M (D) 8-12, 13, 14, 15, 16
 W (EE) 7-12, 13, 14, 15, 16
 4W (4E) 7-12, 13, 14, 15, 16
 6E 7-12, 13, 14, 15, 16



Walker II

40784-14 **Black Pebbled Leather**
 40784-64 **Brown Pebbled Leather**
 N (B) 9-14
 M (D) 8-14, 15, 16
 W (EE) 8-14, 15, 16
 4W (4E) 8-14, 15, 16
 6E 8-14, 15, 16

HCPCS Code A5500 Traveler, Traveler V, Victor, Walker II, Journey II

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable, Moldable, Dual-density Insert
- Extended Medial Heel Stabilizer
- Lightweight, Flared Polyurethane Outsole
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Rocker Bottom with Wide Shank
- Tempered Steel Shank
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

66



Doubler

40822-11 **Black Pebble Leather**

HCPCS Code A5500 Doubler

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Moldable, Total Contact, Dual-density Insert
- Seamless, Moldable Deerskin Lining
- Polyurethane Outsole with Wide Shank
- Tempered Steel Shank
- Extended Medial Heel Stabilizer
- Accommodates Prescribed Orthotics

N (B) 9-13,14
 M (D) 7-14, 15, 16
 W (EE) 7-14, 15
 4W (4E) 7-14

Drew Shoe uses genuine Velcro®.

Men's

Men's Athletic



Force

40960-14 Black Leather

40960-22 White Leather



Force V

44714-14 Black Leather

44714-22 White Leather



Surge

40760-19 Black Leather/Nubuck/Mesh

40760-21 White Leather/Nubuck/Mesh

40760-14 Black Leather

40760-43 Grey Leather/Nubuck/Mesh

HCPCS Code A5500 Force, Force V, Surge, Pulse



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Tri-support Insole
- Rocker Bottom and Wide Shank with Built In Heel Cushion Pad
- Durable Rubber Outsole and EVA Midsole
- Tempered Steel Shank
- Extended Medial Heel Stabilizer
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics

*N (B) 9-13, 14 *No narrow width for Pulse

M (D) 7-13, 14, 15, 16

W (EE) 7-13, 14, 15, 16

4W (4E) 7-13, 14, 15, 16

6E 7-13, 14, 15, 16

401



Pulse

40794-14 Black Calf

Men's

Men's Athletic



Voyager

40890-14 Black Leather (color shown)
40890-22 White Leather

**NEW
Color!**



Aaron

40893-02 White Leather/Mesh with Grey
40893-04 White Leather/Mesh with Blue
40893-09 Black Leather/Mesh with Silver
40893-19 Black Leather/ Black Nylon Mesh New Color! (color shown)
40893-3V Cream Suede/Mesh
40893-43 Grey Suede/Mesh
40893-4V Navy Suede/ White Mesh (color shown)
40893-68 Dark Brown Leather/ Nylon

HCPCS Code A5500 Aaron

DD PF RI



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Drilex® Lining with AEGIS Microbe Shield™ for Odor and Stain Control
- Flared, Lightweight Rubber Outsole with EVA Wedge Midsole
- Removable, Polyurethane Footbed with Heel Cushion and Arch Support
- Firm Heel Counter
- Full Grain Leathers
- Wide, Synthetic Fiber Shank
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics



N (B) 9-13, 14
M (D) 7-13, 14, 15
W (EE) 7-13, 14, 15
4W (4E) 7-13, 14, 15
6E 7-13, 14, 15



Jimmy

44935-14 Black Leather
44935-22 White Leather
44935-32 Bone Leather (color shown)



Lightning II

40805-07 White Leather / Grey Mesh
40805-19 Black Leather/ Black Mesh
40805-21 White Leather/ White Mesh (color shown)

HCPCS Code A5500 Voyager, Jimmy, Lightning II, Lightning II V

DD PF RI

- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- Removable, Polyurethane, Drilex® Lined Insole
- Breathable Soft Fabric Lining
- Durable Rubber Outsole and EVA Midsole
- Dual Shank System
- Wide, Synthetic Fiber Shank
- External TPU Shank
- Firm Heel Counter
- Rocker Bottom
- Soft, Foam-padded Tongue and Collar
- Accommodates Prescribed Orthotics



N (B) 9-13, 14
M (D) 7-13, 14, 15, 16
W (EE) 7-13, 14, 15, 16
4W (4E) 7-13, 14, 15, 16
6E 7-13, 14, 15, 16




Lightning II V

44735-07 White Leather / Grey Mesh (color shown)
44735-19 Black Leather/ Black Mesh
44735-21 White Leather/ White Mesh

Men's

Drew Shoe uses genuine Velcro®.
Velcro® is a registered trademark of Velcro Industries B.V.

THERMOLITE®



Rockford

- 40808-18 Black Tumbled Leather
- 40808-68 Brown Tumbled Leather
- 40808-9R Wheat Nubuck

THERMOLITE®



Boulder

- 40920-18 Black Tumbled Leather
- 40920-68 Brown Tumbled Leather

Trail

- 40866-2U Olive Suede Leather
- 40866-4V Navy Suede Leather



HPCPS Code A5500 Boulder, Rockford, Trek, Trail

DD PF RI SS



- Plus Fitting System® - Two Removable Footbeds for Added and Double Depth
- **Slip-resistant, Oil Resistant, Non-marking Rubber Outsole**
- Waterproof Leather Upper
- Tempered Steel Shank
- Accommodates Prescribed Orthotics

Rockford/Boulder Sizes

- M (D) 8-13, 14, 15, 16
- W (EE) 8-13, 14, 15, 16
- 4W (4E) 8-13, 14, 15, 16
- 6E 8-13, 14, 15, 16



Trail/Trek Sizes

- N(B) 9-13
- M (D) 8-13, 14, 15, 16
- W (EE) 8-13, 14, 15, 16
- 4W (4E) 7-13, 14, 15, 16
- 6E 7-13, 14, 15, 16



Trek

- 40697-15 Black Nubuck
- 40697-68 Dark Brown Leather

